



Determination of contaminating particles in cell therapy products; a proof of principle study.

Paul Dyer, Karessa White, Dikran Khachadourian, Bernardo Cordovez
Halo Labs, Burlingame, CA (paul.dyer@halolabs.com)



Introduction

Distinguishing between different subvisible particles (SVP) species within a cell therapy product is not facile and has proven challenging using standard techniques such as flow cytometry^[1]. This is, in part, due to the cell therapy product contributing to the SVP population.

This has recently been exemplified by the recent FDA issuance of a Form 483 concerning particle contamination of the Car-T cell therapy Kymriah® with wood, cellulose, brass, and steel^[2]. It is generally considered to be attributed to the cryopreservation bags. This recent observation has highlighted a need for a new approach to determining SVPs within cell therapies.

Using a novel technique, we employed backgrounded membrane imaging (BMI), side illumination membrane imaging (SIMI), and fluorescence membrane microscopy (FMM), available on the Aura+, to examine the presence of plastic and other extrinsic particles in a model Car-T cell population.

This proof of principle study aimed to determine the viability of extrinsic particle analysis using Aura+ in complex cell therapies.

BMI and SIMI Explained

Briefly, subvisible particles were immobilized on a backgrounded 0.8µm porous filter membrane using vacuum (15 inHg) and imaged using epi-illumination (BMI) and side-angle illumination (SIMI) to determine all particles and extrinsic particles, respectively. Extrinsic particles remain rigid and thus exhibit a positive SIMI signature (see figure 1).





Figure 1. Schematic representation of BMI and SIMI.

a) BMI utilizes epi-illumination to capture all particles in high contrast; b) SIMI can distinguish rigid extrinsic particles from flat biologic particles.


Extrinsic Particle Counting

To model the presence of extrinsic SVP in cell therapies, Car T cells were mixed with the NIST particle standard (abraded Ethylene tetrafluoroethylene)^[3], cellulose fibers suspended in filtered water for injection (WFI), and metallic particles (commercially available iron fillings suspended in filtered WFI). To assess particle populations, 50µL of each sample was loaded in triplicate onto a backgrounded filter membrane (0.8 µm pore size) and analyzed. Herein, a mean extrinsic particle count of $1.1 \times 10^4/\text{mL}$ (SD±0.05), $2.9 \times 10^4/\text{mL}$ (SD SD±0.15), and $0.5 \times 10^4/\text{mL}$ (SD±0.05) respectively (see figure 2).



Residual Dynabeads™ Detection


Identifying residual Dynabeads* in cell therapies, less than 100 per 3×10^6 cells, is an important consideration^[4]. SIMI analysis can be used to count residual beads in a Car T cell preparations. Herein, we analyzed $1 \times 10^6/\text{mL}$ Car-T cells mixed with CD3/CD28 activation Dynabeads (figure 3).



*Disclosure: Dynabeads refers to the magnetic beads produced by Thermo Fisher Scientific, Inc. Halo Labs is not affiliated with Thermo Fisher Scientific, Inc., and references to Dynabeads or any other third-party trademark do not imply sponsorship, endorsement, or approval.

FMM for multiplex analysis

The Aura+ incorporates three fluorescence channels, FL-1 (ex:440 em:500), FL-2 (ex:376 em:440), FL-3 (ex:482 em:524). Combined with BMI and SIMI analysis, FMM allows for total particle analysis within a complex particle cell mixture. Briefly, Car T cells ($5 \times 10^5/\text{mL}$) were treated with 5x freeze-thaw cycles (to stimulate cell death), stained with DAPI ($2\mu\text{g}/\text{mL}$) to resolve nuclear details, and Thioflavin T (ThT) to proteinaceous particles. Four particle subtypes were identified; plastic (SIMI⁺), non-viable cells (Dapi^{high}), viable cells (Dapi^{low}), and protein aggregates (ThT⁺ Dapi⁻) (figure 4).



Conclusion

This proof of principle study highlighted the capability of Aura+ to identify, count, and size different subvisible particle populations within a complex mixture of cells, protein, and plastic particles.

Thus, providing an important breakthrough in counting contaminating extrinsic particle populations within cell therapy products.

References

- [1] Molina *et al.* Cytotherapy 24 (12) Dec 24, (2022) DOI: 10.1016/j.jcyt.2022.07.014
- [2] <https://www.fda.gov/media/174213/download>
- [3] Ripple *et al.* NIST SP (2019) DOI: 10.6028/NIST.SP.260-193
- [4] White *et al.* Journal of American college of Toxicology 14 (4) (1995) DOI: 10.3109/10915819509008701
- [5] Foglieni *et al.* Histochem Cell Biol 115 (2001). DOI: 10.1007/s004180100249